

# ALMICRO

## BOD Incubator



### **Specification:**

Double wall construction. Inner chamber made of anodized aluminum or highly polished stainless steel sheet. Supplied with two/ three removable perforate shelves. Exterior made out of thick mild steel duly finished in white stoving enamel/ powder coated paint with mat finished color combinations. Outer double wall door is provided with magnetic gasket and lock & key arrangement. Inner door unbreakable acrylic transparent sheet for inspection of material in the chamber without disturbing the chamber temperature which is controlled by Deluxe Digital Controller from 5°C to 50°C 0.5°C. Standard cooling unit with a sealed compressor with protective devices has been provided below the chamber for safe & efficient working. Heating arrangement has been made within the chamber for higher temperature working. Inner air circulating system is provided to achieve uniform temperature. A control panel is provided with various indicators, on – off switch. Supplied with cord and plug etc. to work on 220V AC 50Hz single Phase.