

# ALMICRO

## Bacteriological Incubator



### **ALMICRO BACTERIOLOGICAL INCUBATOR:**

Standard double wall fabrication. Inner chamber made out of Richy anodized aluminum or highly polished stainless steel sheet and outer made out of thick mild steel sheet finished in white staving enamel paint or mat finished pastel color combinations. The gap is filled with special glass wool for temp. Insulation. Double wall door with double glass window for observation in the chamber.

Temperature range  $5^{\circ}\text{C}$  above ambient to  $80^{\circ}\text{C} \pm 2^{\circ}\text{C}$  controlled by a thermostat L shaped thermometer is provided for taking temp. readings. Panel is provided with pilot lamps for mains & heaters. On off switch temp. control knob etc. The equipment is workable on 220 V AC 50Hz single phase. Heating is achieved from 3 sides.

Perforated removable shelves also provided.