

ALMICRO

Hot Air Oven



Standard double wall fabrication. Inner chamber made of richly anodized fabricated out of thick mild steel sheet duly finished in stored enamel off-white paint or pastel mat finish colors combinations. The gap between the walls is filled with special grade grass wool for heat insulation. Inner chamber is ribbed for placing shelves at convenient levels. Provided with two or three remove able shelves. Special elegant door closing device has been provided. Nichrome wire heating elements provided on three sides attain quick and uniform heating in range of 50°C to $250^{\circ} \pm 1^{\circ}\text{C}$ controlled by Microprocessor based PID Digital Temperature controller fitted on front door) Air circulation fan is provided for uniform heating. Supplied with cord and plug. The equipment is suitable to operate on 220V AC 50 Hz single phase.