



ALMICRO Inclined Metallurgical Microscope comes fitted with Monocular observation tube inclined at 45°. It is provided with Bright field horizontal illuminator with centering provision, field and aperture diaphragm, and slots for dropping filter. It comes equipped with Quadruple revolving nose piece. There is a Special graduated mechanical stage for easy study of the samples and is supplied with variable light control transformer and following optics in box :

**LWD Objectives**  
**Wide Field Eye pieces**

M10x and M45x (SL)  
WF10x and H15x



**NSIC**

