

ALMICRO[®]

Student Polarizing Microscope PM-9



Features :

Limb type hinged body inclinable 90°.

Bertrand lens, sliding analyser 1/4 wave & full wave compensator.

Ball Bearing mounted, rotatable stage divided 0-360° with vernier.


Adjustable condenser N.A. 1.25.

Polariser duly graduated with in-out system.

Objectives P5x, P10x, P40x/45x.

Eye-piece 5x(cross), 10x(ocular) & 15x(optional).

Plano-concave reflector replaceable by sub-stage L.E.D. Lamp(optional)



NSIC

