

ALMICRO

SEROLOGICAL WATER BATH



Specification:

Double wall construction inner being made out of highly polished thick stainless steel and outer made out of thick mild steel sheet finished in staving white enamel paint or staple mat color combinations. The gap is filled with special glass wool for temp. insulation. The bath is provided with a drain cock to easily drain out the water and clean the chamber when required. A pyramidal shaped cover is provided.

Immersion type heating elements are fitted at bottom. Temp. range 5^oC above ambient to 80^oC with accuracy $\pm 0.5^{\circ}\text{C}$ to 1^oC is controlled by capillary thermostat regulated by a graduated knob on the panel, which is also provided with pilot lamp. Etc. The equipment is workable on 220V AC 50 Hz single phase. The baths are supplied without racks. Thermometer and stemming arrangement.