


	BM-3 ISI Student Microscope
 Eye-Piece	Wide-Field 10X Eye-piece or set of 10X & 15X Huygenian Eye-pieces.
Nosepiece	Triplet revolving nosepiece with click stops
Stage	High Quality S.S. Stage clips provided to hold the slide.
Focusing	By vertical movement of the observation tube, with coarse adjustment using rack-pinion movement and fine motion by screw level mechanism.
Condensor	Abbe condenser having Numerical aperture of 1.25 with iris diaphragm and swing out filter holder. This condenser unit incorporates high efficiency optical system for optimum utilisation of light.
Illuminator	Arc-Reflector set for using natural light to study the slide.
Objective	Achromatic 10X and 45X.
Packing	In mould styrofoam box or wooden storing cabinet with lock and key (on demand)
Accessories	Vinyl dust cover, Operating manual



Model : BM3 - ISI

ALMICRO[®] Compound Student Microscopes have well designed heavy and sturdy pressure die casted monocular body inclinable to 90° to meet the requirements of School and college students