

ALMICRO[®]

Binocular Polarizing Microscope PM-11



Features :

- Binocular head replaceable with Monocular, and Trinocular head.
- Bertrand lens, Sliding Analyzer 1/2 & 1/4 wave compensator.
- Rotatable stage divided into 0-360° with vernier.
- Adjustable Condensor N.A. 1.25.
- Polarizer duly graduated with in-out system.
- Built in light arrangement with 6V20W halogen lamp.
- Solid State light control arrangement.
- Paired eye piece 5x(cross), 10x(ocular), 15x (optional).
- Objectives P5x, P10x, P40x/45x



NSIC

