

ALMICRO

Water Bath (12 Holes)



Standard double wall construction, inner being of stainless steel outer of m.s. duly staving paint finish and gap filled for temp. insulation with glass wool. The cover with 75mm (3") holes with concentric rings. Temp. range 5⁰C above ambient to 95⁰C controlled by a thermostat. Heating is attained by an immersion type element with accuracy of $\pm 1^{\circ}\text{C}$.

Supplied complete with plug & cord. Panel provided with pilot lamps, thermostat control knob etc.

