

Pathological Monocular Microscope  
ISI Marked to IS:4381/1967  
equipped with:

- ▶ One 10X Wide Field Eye -piece.
- ▶ 4 nos. Objectives, Parfocal Achromatic, 4x, 10x, 40x SL and 100x SL Oil immersion.
- ▶ Stage : Horizontal Mechanical of size 120 X 125mm min.
- ▶ Illumination : 6V,20W Halogen lamp with on/off switch, light intensity regulator and Plano concave mirror attachment suitable for 220V, 50 Hz. AC supply.



## Micro Measures & Instruments

#444, Arya Nagar, Jagadhri Road, Ambala Cantt - 133001, India

Ph: +91-9896051072, 9416447169 +91-171-2663258, 2602150, 4005258.

web : [www.almicromicroscopes.com](http://www.almicromicroscopes.com) mail : [sales@almicroinstruments.com](mailto:sales@almicroinstruments.com)