

# ALMICRO<sup>®</sup> Ore Microscope



## Features :

**Both Incident & Transmitted lights are provided for ore Microscopy**

Binocular head replaceable with Monocular, and Trinocular head.

Bertrand lens, Sliding Analyzer 1/2 & 1/4 wave compensator.

Rotatable stage divided into 0-360° with vernier.

Adjustable Condensor N.A. 1.25.

Polarizer duly graduated with in-out system.

Built in light arrangement with 6Y20W halogen lamp.

Solid State light control arrangement.

Paired eye piece 5x(cross), 10x(ocular), 15x (optional).

Objectives P5x, P10x, P40x/45x



**NSIC**

