

# Junior Student Microscope

(Model : BM-2)

Junior Student Microscope  
equipped with:

- Eye Piece** One 10X & one 15X Hygeniun Eye-Piece
- Nosepiece** Triple revolving nosepiece with click stops.
- Stage** High quality S.S. Stage clips provided to hold the slide.
- Focusing** By vertical movement of the Observation tube, with coarse adjustment using rack pinion movement and, fine motion by screw lever mechanism,
- Condenser** Abbe condenser having Numerical aperture of 125 with iris Diaphragm and swing out filter holder. This condenser Unit incorporates high efficiency optical system for optimum utilization of light'
- Light** arc reflector set for using natural light to study the slide
- Objective** Achromatic 10x and 40x.
- Packing** in mould Styrofoam box

