

BM-8bi Pathological Binocular Microscope with co-axial focussing control

Pathological Binocular Microscope
with co-axial focussing control
ISI Marked to IS:4381/1967
equipped with:

- ▶ Two 10X Wide Field Eye -piece.
- ▶ 4 nos. Objectives, Parfocal Achromatic, 4x, 10x, 40x SL and 100x SL Oil immersion.
- ▶ Stage : Horizontal Mechanical of size 120 X 125mm min.
- ▶ Illumination : 6V,20W Halogen lamp with on/off switch, light intensity regulator and Plano concave mirror attachment suitable for 220V, 50 Hz. AC supply.
- ▶ Eye Pieces for Binocular Microscope with coaxial focussing control confirm to IS : 8275:1976



Micro Measures & Instruments

#444, Arya Nagar, Jagadhri Road, Ambala Cantt - 133001, India

Ph: +91-9896051072, 9416447169 +91-171-2663258, 2602150, 4005258.

web : www.almicroscopes.com mail : sales@almicroinstruments.com