

ALMICRO[®] *ultra* Series

BM-3 Student Microscope

◆ Eye-piece	Wide field 10X Eye-piece specially designed for ultra series.
◆ Nosepiece	Triplet revolving nosepiece with click stops.
◆ Stage	High quality S.S. Stage clips provided to hold the slide.
◆ Focusing	By vertical movement of the Observation tube, with coarse adjustment using rack-pinion movement and, fine motion by screw lever mechanism.
◆ Condenser	Abbe condenser having Numerical aperture of 1.25 with iris Diaphragm and swing out filter holder. This condenser unit incorporates high efficiency optical system for optimum utilization of light.
◆ Light Source	Battery operated LED pin-type Lamp with on and off switch and an optional arc-reflector set for using natural light to study the slide.
◆ Objective	Achromatic 10x and 40x.
◆ Packing	In mould Styrofoam box or wooden storing cabinet with lock and key (on demand)
◆ Also Included	Vinyl dust cover, lens cleaning wipes (100 no.s), operating manual.

