

ALMICRO[®] *ultra* Series

BM-4 Junior Medical Microscope

Eye-piece	Wide field 10X Eye-piece specially designed for ultra series.
Nosepiece	Triplet revolving nosepiece with click stops.
Stage	High quality S.S. Stage clips provided to hold the slide.
Focusing	By vertical movement of the Observation tube, with coarse adjustment using rack-pinion movement and, fine motion by screw lever mechanism.
Condenser	Abbe condenser having Numerical aperture of 1.25 with iris Diaphragm and swing out filter holder can be moved through rack-pinion. This condenser unit incorporates high efficiency optical system for optimum utilization of light.
Light Source	Battery operated LED pin-type Lamp with on and off switch and an optional arc-reflector set for using natural light to study the slide.
Objective	Achromatic 10x and 40x.
Packing	In mould Styrofoam box or wooden storing cabinet with lock and key (on demand)
Also Included	Vinyl dust cover, lens cleaning wipes (100 no.s), operating manual.

