

## PORTABLE AUTOCLAVE



Suitable for general laboratory use sterilization. Stainless steel or aluminium vertical body having no joints. Pressure range 15psi to 18psi. Domed cover can be sealed by a jointless gasket which is clamped to the vessel by wing nuts by hinged bolts which engage. Spring loaded safety valve has threaded cap for pressure adjustments. Supplied complete with plug and cord and suitable to work on 220V AC 50 Hz single phase.



**Micro Measures & Instruments**



444, Arya Nagar, Ambala Cantt, Haryana - 133001, INDIA



+91-9896051072, +91-9416447169, +91-171-2663258, +91-171-4005258



sales@almicroinstruments.com, almicroinstruments@gmail.com



www.almicroinstruments.com, www.metallurgicallaboratoryinstruments.com