

ALMICRO AUTOCLAVE



Autoclave Vertical: -

A versatile, electrically heated sterilizer under saturated steam of adjustable pressure of 10 to 20 psi. Double walled construction with boiler and steam jacket made of thick stainless steel and gap is filled with high grade glass wool to bring the temperature loss to the minimal level. The outer chamber is made out of thick mild steel. The lid is made of thick stainless steel plate (joint less) or gunmetal plate and can be sealed by a joint less neoprene rubber gasket by wing nuts through hinged bolts by engaging slotted lugs or rim of the cover.

Provided with dial pressure gauge, air/steam release cock, spring loaded safety valve adjustable to any required point from 10 psi to 20 psi with accuracy ± 3 psi. A drain cock is also provided. A double safety valve is provided along with water inlet. Supplied complete but without dressing bins.