


ADVANCE ROTARY MICROTOME

(Manual microtome)



SALIENT FEATURE

- SECTION THICKNESS DISPLAY
- ADJUSTABLE KNIFE ANGLE
- COARSE FREE HAND-WHEEL
- SECTION THICKNESS ADJUSTABLE KNOB
- SPACIOUS SECTION WASTE TRAY
- INTEGRATE FINGER GUARD
- CROSS ROLLER GUIDE WAYS



PRODUCT DESCRIPTION:

Advance rotary (manual) microtome is designed for manual routine paraffin sectioning but can also be used for cutting harder material. It is perfectly suited for those who expect precision sectioning and appreciate the feel and extra measure of control during the feel and extra measure of control during manual sectioning.

TECHNICAL SPECIFICATION:

- Section thickness range: 0-60 μ m
- Setting values: From 0-5 μ m-2 μ m in 0.5 μ m increments
From 2 μ m-20 μ m in 1 μ m increments
From 20 μ m-60 μ m in 5 μ m increments
- Horizontal specimen stroke: 28mm
- Vertical specimen stroke: 60mm
- Precision error: \pm 5%
- Sample Orientation: 8° along the X-Y axis, routable 360°
- Maximum section size: 50x45mm
- Dimension incl. tray: 300x570x270mm