



ALMICRO DS-50 Stereo Microscope comes equipped with inclined stereo body, Binocular observation Head having adjustable interpupillary Distance, and diopter correction eye piece tube.

This Microscope is provided with built in illumination both for Incident (EPI) light and Substage (transmitted) light for individual or simultaneous use.

Suitable transformer is fitted in base for 6v-15w lamps.

It has a Revolving Turret Mount for Quick Change for objective lenses.

It is Supplied with following optics :

Eye piece - WF 10x (Paired), Optional -15x Eye Piece

Objectives - 2x and 4x

 **Micro Measures & Instruments**