

# ALAMICRO HDMI CAMERA

## Plug and Play Microscope CMOS USB 16MP Camera



### Specification

- 16 megapixel CMOS image sensor;
- 1/2.33 inch image sensor;
- HDMI 1080P 60 frame output;
- USB2.0 1080P 30 frame output;
- 4K video recording
- 16 megapixel photos, photo size 4608x3456;
- Freezing/defrosting function;
- Electronic zoom;
- Support manual/automatic exposure/one-button exposure, exposure compensation adjustment;
- Supports manual/automatic white balance/one-button white balance, red, green and blue adjustable;
- The function of 8 drawing lines is provided. The line color, thickness and position can be adjusted;
- center cross display/close;
- Support TF card storage, maximum support of 64GB;
- The remote control operation, the camera does not shake, the image is clearer;

