



ALMICRO Projection Microscope (KW-500)

This Microscope is meant for group discussion & teaching a group of students.

SPECIFICATIONS :

- ILLUMINATION : It comes equipped with 12V 100W Halogen Lamp.
- STAGE : 110 x 110mm square stage is provided with stage clips to hold the slide (optional Mechanical Stage).
- FOCUSING : Focussing control knob on both sides of the instrument is conveniently located for easy operation of the rack and pinion with 15mm range.
- NOSE PIECE : Triple revolving nosepiece with positive click stops.
- SCREEN & WALL PROJECTION : 160mm dia. optically true ground glass screen is provided to have high resolving power. The screen housing is hinged to the body frame and can be removed to one side to have projection on the wall or screen in a dark room.
- OBJECTIVES : This Microscope is equipped with three Achromatic objectives 5x, 10x & 20x.
- PACKING : This instrument is supplied in sturdy cardboard box with styrofoam support.

Micro Measures & Instruments

