



ALMICRO Student Projection Microscope (Doom Type) (KW-700C)

This Microscope is Compact in Size & can easily project distortion free clear image on a fine ground glass screen.

SPECIFICATIONS :

- Body : Monocular fixed straight tube (Length - 160mm)
- Nosepiece : Triple revolving Nosepiece with positive click stops
- Stage : 110 x 110mm square stage is provided with stage clips to hold the slide (optional Mechanical Stage).
- Illumination : LED Lamp
- Objectives : Achromatic 4x, 10x & 40x
- Eyepiece : Hyugenanian 10x
- Screen : 150mm diameter glass projection screen is provided which can be rotated over 360 degrees on Microscope tube.
- Packing : Securely packed in Cardboard & Styrofoam box.

Micro Measures & Instruments

