

**ALMICRO**

## **DIGITAL CENTRIFUGE MACHINE (KW-80A)**



### **Applications:-**

- The best model for clinical, Pathological, Biochemistry and pharmacy/pharmaceutical Labs and to separate blood serum

### **Specifications:-**

- General Purpose Bench Top laboratory Centrifuge Machine.
- Heavy duty C. R. C construction finished with powder coating.
- See through Acrylic moulded lid which is light weight to avoid accident.
- Inter lid lock and metallic centrifuge bowl for safety.
- Microprocessor based programmable speed control from 500 to 5000 r.p.m depending upon the type of rotor and Digital timer (0-59 minutes)
- Supply is completed with 8x15ml Angle head, instruction manual, spare fuses, cord & plug to work on 220volts 50Hz A.C.

## **Micro Measures & Instruments**

#444, Arya Nagar, Ambala cantt - 133001  
Customer care : +91 8929310310

Email : sales@almicroinstruments.com  
www.almicroinstruments.com