

# TOOLMAKERS MICROSCOPE (MODEL : KW-900)



This Microscope has X-Y movement on ball bearing slides, having micrometer (50mm dia) Least Count - 0.01 or 0.005mm with goniometer. Fixed magnification 30x. 150x150mm stage, and 360° rotary head. Having both transmitted & incident light. This comes packed in a wooden box.

## Technical Specifications:

- 🔍 **Observation Tube**  
Monocular, Inclined at 30-degrees.
- 🔍 **Base**  
Large and heavy base provides extra overall rigidity to the instrument.
- 🔍 **Stage**  
150mm X 150mm size assembled on ball bearing guides to provide accurate and smooth stage stroke upto 50mm in each direction.
- 🔍 **Illumination**  
Sub-stage lamp emits transmitted light from a bottom source equipped with halogen lamp.

## 🔍 Micro Measures & Instruments

📍 444, Arya Nagar, Ambala Cantt, Haryana - 133001, INDIA

☎ +91-9896051072, +91-9416447169, +91-171-2663258, +91-171-4005258

✉ sales@almicroinstruments.com, almicroinstruments@gmail.com

🌐 www.almicroinstruments.com, www.metallurgicallaboratoryinstruments.com