

X-555 Plus With the features of Infinite Optical System, Effective Illumination LED Pointer and Images Coherence, It is Widely used In Scientific Research and Teaching Demonstration.



X-555 Plus

Optical System	Infinite Optical System
Viewing Head	Compensation Free Trinocular Head, Inclined at 30°, 360° Rotatable. Interpupillary Distance 48-75mm (1PC) Compensation Free Binocular Head, Inclined at 30°, 360° Rotatable. Interpupillary Distance 48-75mm (4
Eyeiece	Extra Wide Field Eyepiece EW 10x80 (10PCS)
Nosepiece	inverse Quintuple Nosepiece
Objective	Infinite Plan Achromatic Objective 4x, 10x, 40x or 100x
Condenser	Swing Condenser N.A.0.9/0.25
Focusing	Coaxial Coarse & Fine Focus Adjustment, Fine Division 0.001mm
Stage	Double Layers Mechanical Stage 185x142mm, Moving Illumination, Halogen Lamp 24V/100W
Kohler Illumination Pointer	External Illumination, Halogen Lamp 24V/100W Green LED Pointer, Brightness Adjustable

Micro Measures & Instruments