

RESEARCH
MICROSCOPES
x-series

PENTA HEAD MICROSCOPE



Technical Specifications

Viewing Head	Sixing Binocular Head, Inclined at 45° (5pcs)
Eyepiece	Wide Field Eyepiece WF 10x (10/vox)
Nosepiece	Quadruple Nosepiece
Objective	Achromatic Objective 4x, 10x, 40x or 100x
Condenser	Abbe Condenser NA 1.25
Focusing	Coaxial Coarse & Fine Focus Adjustment, Fine Division 0.002mm
Stage	Double Layers Mechanical Stage 140x140mm/75x50mm
Illumination	S-LED Illumination, Brightness Adjustable Halogen Lamp 6V/20W Halogen Lamp 42V/100W
Pointer	Green LED Pointer, brightness Adjustable

X 555 Penta head Microscope is mainly used for teaching A demonstration purpose They can also be used for biological analysis for more experts at the same time. They can be widely used in medical and sanitary established laboratories, institutes, colleges and universities.