



SLIDE STAINING MACHINE



SALIENT FEATURE

- FEATHER TOUCH KEYPAD
- DIGITAL DISPLAY
- BATTERY BACKUP TO RETAIN DATA
- EASY OPERATED
- COMPACT SIZE

PRODUCR DESCRIPTION:

SLIDE STAINING MACHINE operated through a keyboard. The most advance Slide staining machine to enhance control of different requirement in staining technique starting from the routine haematological techniques to the most advance of staining in cytology and histology. It cut down the time and labour required to manually dip the slides through different reagents.

TECHNICAL SPECIFICATIONS:

- ✓ Digital Display for adjusting the timing sequence
- ✓ Battery backup to retain data
- ✓ Durable construction
- ✓ Twelve staining troughs
- ✓ High Specimen Throughput
- ✓ Easy operated and Easy Maintenance
- ✓ Compact size
- ✓ Maximum output by offering high specimen
- ✓ Minimum reagent exposure
- ✓ Easy Maintenance
- ✓ Trough Volume- 450ml
- ✓ Trough Capacity- 15 slides
- ✓ Ideal for haematology and staining with 12 Glass staining trough, 1 slide carrier
Dust cover, Oil can, Manual Construction
- ✓ Suitable for work on 220V, Single phase, 50Hz, AC supply

INSTRUMENTS WILL BE SUPPLIED IN BOX WITH STYROFOAM PACKING.