



Spectro Polishing Machine

Avail from us an exclusive range of quality Spectro Polishing Machine that is made out of latest technology and quality raw material procured from reliable sources. This machine is compact in size & robust in design. Paper replacement is easy. Built in dust collector, collects all the hazardous particles in a separate disposable bag. We meet customize demands at really affordable prices

FEATURES

DISC COVER

Machine to Suit the disc housing & to cover the disc plate by 180° for safety

DISC PLATE

Made up of Steel with air gripping arrangement to aluminum disc and with locking system to keep the aluminum disc in position

DUST COLLECTION UNIT

At the back of the disc housing the pipeline connected to the dust collecting unit

NOISELESS OPERATION

In continuous process, there is noise less operation

PROTECTION

Anti fungal, rust proof paint used for Extra protection

SPECIFICATIONS

MAIN MOTOR

1hp, 1400rpm

POWER SUPPLY

440/460 volts, 3 phase, AC Supply

AIR GRIPPING MOTOR

2HP, 2800rpm, L

STARTERS

L and T make push button type for main and air gripping motor