


TISSUE EMBEDDING STATION



SALIENT FEATURES

- LED DISPLAY
- EASE OF OPERATION
- HIGH STANDARD OF CONVENIENCE
- STAINLESS STEEL
- COMPACT DIMENSION



TECHNICAL SPECIFICATION:

- Advanced flexible heating method, to overcome the traditional shortcomings of excessive heating temperature, quick heating and high precision, as well as double overheat protection, safe, reliable and energy saving.
- Tissue Embedding Station characterized by its Compact Dimension, Ease of Operation and High Standard of Convenience.
- Large LED display allows user to see temperature without scrolling.
- Left and Right thermal Chambers, Hot plate, Forceps wells, Dispenser all individually temperature controlled.
- Light indicators keys
- Heated Holder for 6 Forceps
- Stainless Steel Paraffin large tank of 3-5 liters capacity