

ALMICRO VDRL SHAKER



ALMICRO ROTARY SHAKER AM-131

(VDRL Rotator) Fixed Speed

Suitable for VDRL and other agglutination tests, blood grouping tests and for mixing of solutions in small flasks, beakers or bottles. Body made of sturdy mild steel sheet duly finished in stoving enamel paint. Tray type platform size 300 x 300mm can take cavity slides. Bottles flasks etc, for shaking at a fixed speed by the help of spring wires. Standard fixed speed 180 RPM. An automatic 60 minutes timer can be provided against extra cost as optional by which shaking duration can be set. To work on 220V AC 50 Hz single phase.

NOTE: Variable speed rotator can be offered on demand.



MICRO MEASURES & INSTRUMENTS



+91-9896051072, +91-9416447169, +91-171-4005258



sales@almicroinstruments.com, almicroinstruments@gmail.com



www.almicroinstruments.com, www.metallurgicallaboratoryinstruments.com