

TRAVELLING MICROSCOPE (MODEL : VM-1)



Our Travelling microscope (**VM-1**) consists of a cast iron base with machined vee-top surface and is fitted with three levelling screws, one being fixed the microscopes can travel Horizontally 22 cm. and vertically 15 cms. with the help of slides. Two slow motion knobs are provided for taking accurate readings. Vernier readings 0.01 mm. or 0.02 mm. Milky cenolite sheet is provided on the base for the placing objects on horizontal stage. Nicely packed in wooden box.

with following optics:

- Objective – Achromatic 2x (FL-50mm)
- Eye Piece – 10x with cross hair line

www.almicroinstruments.com