

ALMICRO[®]

DELUXE
TRAVELLING MICROSCOPE
(MODEL : VM-2B)



Deluxe Travelling Microscope (VM-2B) for **three axes measurement** Horizontal, vertical and right angled transverse motion with rack and pinion mechanism. The microscope can travel horizontally 20cm, vertically 15cm & right angled transverse movement 25mm. Nicely packed in wooden box.

with following optics:

- Objective – Achromatic 2x (FL-50mm)
- Eye Piece – 10x with cross hair line

www.almicroinstruments.com

ALMICRO[®]
MICROSCOPES



NSIC

