


FULLY AUTOMATIC MICROTOME



SALIENT FEATURES

- ERGONOMIC DESIGN
- LCD DISPLAY
- SOFT TOUCH
- WASTE TRAY FOR EASY CLEANING
- MULTIPLE MODES FOR SECTIONING



PRODUCT DESCRIPTION:

The latest technology and perfectly harmonizes with an ergonomic design and user safety.

TECHNICAL SPECIFICATION:

- ❖ Rust Proof Construction
- ❖ Lockable hand wheel at any position allowing change of blade/ specimen easier and safe
- ❖ Large storage capacity on the top of the instruments for easy access of tools
- ❖ Ergonomic and user friendly design
- ❖ Touch fast forward/ backward function
- ❖ Slicing/ Trimming mode
- ❖ LCD display with feather touch keypad
 - Horizontal Feed: 28mm
 - Vertical Specimen stroke: 70mm
 - Specimen holder rotation at angle within 360°
 - Specimen clamp rotation at any angle within 360°
 - Specimen Orientation XY-8°
 - Maximum Specimen size 70x70mm
 - Dimensions: 600x350x320mm