

Centrifuge (Doctor Model)

Model No : KW-50



Capacity : 4x15ml, 6x15ml, 8x15ml
(Angle Head - Doctor Model)
Mains Input : 220v,50Hz AC
Max. Speed : Upto 3500 RPM
Speed Control : 5 Step Regulator
Motor : Universal (AC/DC Carbon Brush type)

Centrifuge (Handi Type)

Model No : KW-60



Capacity : 4x15ml, 6x15ml, 8x15ml
(Swing Out Head)
Mains Input : 220v,50Hz AC
Max. Speed : Upto 3500 RPM
Speed Control : 5 Step Regulator
Motor : Universal (AC/DC Carbon Brush type)

Centrifuge (Square Type)

Model No : KW-80



Capacity : 4x15ml, 6x15ml, 8x15ml
(Swing Out Head)
Mains Input : 220v,50Hz AC
Max. Speed : Upto 5000 RPM
Speed Control : Variable speed control knob
& Manual Timer is provided.
Motor : Universal (AC/DC Carbon Brush type)

Centrifuge (Square Type - Digital)

Model No : KW-80 A



Capacity : 4x15ml, 6x15ml, 8x15ml
(Swing Out Head)
Mains Input : 220v,50Hz AC
Max. Speed : Upto 6000 RPM
Speed Control : Variable speed control knob
with digital RPM & digital timer
indicator provided on the panel.
Motor : Universal (AC/DC Carbon Brush type)